Validation liquid – 510 nm

For optical verification of CYANSmart & CYANVision
Store at 2-8°C

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Intended use

The validation liquid is an aqueous solution, used to optically validate a semi-automatic biochemistry analyser at a specific wavelength over the clinically relevant measuring range¹.

Reagent composition

Each kit contains 6 levels of validation liquid, required for the optical validation of one wavelength on one single analyser.

Precautions

The included reagents are environmentally friendly, and do not contain harmful ingredients. Avoid contact with eyes, skin and clothing. Wear laboratory gloves when handling the reagent. In case of eye or skin contact, flush eyes with large amounts of water for several minutes or wash the exposed skin area with water. Use Good Laboratory Practices (GLP) when handling this product.

Preparation

This reagent is ready for use and should be used directly from the glass container. Leave the reagent to reach room temperature (18-30°C) before use.

Storage, stability and disposal

Unopened, this product is stable up to the expiration date specified on the label, when stored at the correct temperature and protected from light and sun exposure. Make sure the reagents are never frozen. Storage temperature is 2-8°C. Do not use the product if deterioration or contamination is suspected or beyond the expiration date or open container stability period. Dispose unused or deteriorated product and waste in compliance with local regulations.

Procedure

1. Ensure that your analyzer is in the best possible working condition (up to date maintenance, autozero, pump calibration...)
2. Perform the "preparation of the analyzer" instructions described in the user manual of the specific analyzer. For the CYANVision (CYANSmart) you find more information in chapter 5 - Daily routine (Chapter 4 in CYANSmart)
3. Adjust the instrument to zero with distilled water.
4. If an RBL measurement is requested, measure distilled water. The RBL value should be between -0,1 AU and 0,01 AU. If not, redo the autozero or consult the user manual of the specific analyzer.
5. Measure level 8 (0,2 AU) 10 times of the validation liquid using the LIN510 program, note down the values.
6. Measure level 7 (0,6 AU) 2 times, note down the obtained values
7. Measure level 6 (1,0 AU) 2 times, note down the obtained values
8. Measure level 5 (1,5 AU) 2 times, note down the obtained values
9. Measure level 4 (2,0 AU) 2 times, note down the obtained values
10. Measure level 3 (2,5 AU) 2 times, note down the obtained values
11. Perform the "end of day" maintenance actions described in the user manual of your analyzer.

Data processing

1. Go to <https://diagnostics.be/check-your-analyzer>
2. Follow the instructions on our website for certification.

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**CYPRESS
DIAGNOSTICS**

REGISTRATION FORM

Laboratory	Name	<i>(Company Name)</i>	
	Coordinates	<i>(Country, City) (E-mail) (Website) (Tel.)</i>	
Analyzer	Execution date		
	Analyzer type	<input type="checkbox"/> CYANSmart (CY009) <input type="checkbox"/> CYANVision (CY014) <input type="checkbox"/> Other: (...)	
	Serial number		

RESULTS

Execution	Level	Measurement
1	Level 8 (±0,2 AU)	<i>(Note result)</i>
2		<i>(Note result)</i>
3		<i>(Note result)</i>
4		<i>(Note result)</i>
5		<i>(Note result)</i>
6		<i>(Note result)</i>
7		<i>(Note result)</i>
8		<i>(Note result)</i>
9		<i>(Note result)</i>
10		<i>(Note result)</i>
1	Level 7 (±0,6 AU)	<i>(Note result)</i>
2		<i>(Note result)</i>
1	Level 6 (±1,0 AU)	<i>(Note result)</i>
2		<i>(Note result)</i>
1	Level 5 (±1,5 AU)	<i>(Note result)</i>
2		<i>(Note result)</i>
1	Level 4 (±2,0 AU)	<i>(Note result)</i>
2		<i>(Note result)</i>
1	Level 3 (±2,5 AU)	<i>(Note result)</i>
2		<i>(Note result)</i>

COMMENTS