

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CF-T1.0AL250V	CL001	Power	Fuse T1,0AL250V 5x20mm (10 pcs)		For CL001-220 On mainboard	N		1:10
CF-T0.5AL250V	CL002, CL007, CL007N, CL014, CM001 & CM001N	Power	Fuse T0,5AL250V 5x20mm (10 pcs)		At the backside (CL002, CL007, CL007N) On mainboard (CL014) For CM001(N) until SN 103592	N		1:10
CL007-S01	CL007 & CL007N	Main Assembly	Rubber ring		Replace when usage is visible	(Y)		
CL007-S02-110	CL007 & CL007N	Main Assembly	Motor - 110V			N		
CL007-S02-220	CL007 & CL007N	Main Assembly	Motor - 220V	<a href="#">Picture / Photo</a>		N		
CL007-S03	CL007 & CL007N	Main Assembly	Driven pulley	<a href="#">Picture / Photo</a>		N		
CL007-S04	CL007 & CL007N	Main Assembly	Plastic gear	<a href="#">Picture / Photo</a>		N		
CL007-S05	CL007 & CL007N	Main Assembly	Plastic fixer for driven pulley			N		
CL007-S06	CL007 & CL007N	Main Assembly	ON/OFF switch			N		
CL007-S07	CL007 & CL007N	Main Assembly	Rubber foot/support		4 per instrument	N		
CL008-S01A-110	CL008	Main Assembly	Motor driver PC board - old, 110V		Until SN 1406..... Not compatible with newer models	N		
CL008-S01A-220	CL008	Main Assembly	Motor driver PC board - old, 220V		Until SN 1406..... Not compatible with newer models	N		
CL008-S01B-110	CL008	Main Assembly	Motor driver PC board - new, 110V		From SN 1409.....on Not compatible with older models	N		
CL008-S01B-220	CL008	Main Assembly	Motor driver PC board - new, 220V		From SN 1409.....on Not compatible with older models	N		
CL008-S02A-110	CL008	Main Assembly	DISPLAY - RPM/TIME PC board (old, 2		Until SN 1101..... Not compatible with newer models	N		
CL008-S02A-220	CL008	Main Assembly	DISPLAY - RPM/TIME PC board (old, 2		Until SN 1101..... Not compatible with newer models	N		
CL008-S02B-110	CL008	Main Assembly	DISPLAY - RPM/TIME PC board (new, 1		From SN 1102000 on Not compatible with older models or CL008N	N		
CL008-S02B-220	CL008	Main Assembly	DISPLAY - RPM/TIME PC board (new, 1		From SN 1102000 on Not compatible with older models or CL008N	N		
CL008-S04A	CL008	Main Assembly	Rubber gasket (old)		Replace when usage is visible Until SN 1308..... Not compatible with newer models Includes height adjustment piece	Y		
CL008-S05A-110	CL008	Main Assembly	Motor- old, 110V		Until SN 1406..... NOT compatible with newer models	N		
CL008-S05B-110	CL008	Main Assembly	Motor - new, 110V		From S/N 1503.... On also compatible with older SN	N		
CL008-S05B-220	CL008	Main Assembly	Motor - new, 220V	<a href="#">Picture / Photo</a>	From S/N 1503.... On also compatible with older SN	N		

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CL008-S05B-220-CK	CL008	Main Assembly	Compatibility kit for new motor for cei		S/N 14090009 - 14090058 (available on request with order of CL008-S05B and mentioning of SN)	N		
CL008-S11A	CL008	Main Assembly	Variable resistance for speed (OLD)		Until S/N 1406.... Not compatible with newer models	N		
CL008-S11B	CL008	Main Assembly	Variable resistance for speed (NEW)		From S/N 1409.... On Not compatible with older models	N		
CL008-S14	CL008	Main Assembly	Transformer for fan			N		
CL008-S18	CL008	Main Assembly	Front label		Do NOT reattach the old label with (super)glue. Reattach with doublesided tape or use a new front label.	N		
CF-T3.15AL250V	CL008 & CL008N	Power	Fuse T3.15AL 250V 5x20mm (10 pcs)		At the back of CL008-220 and CL008N-220	N		1:10
CL008/9-S08	CL008 & CL009	Main Assembly	Power switch		Not compatible with CL008N or CL009N	N		
CL008/9-C-A	CL008 & CL009	Accessories & Consumables	Carbon brush (pair) - old	<a href="#">Picture / Photo</a>	Replace every 500 working hours, Replace every 200h if use is >5-6h/days CL008 until SN 1406... For CL008-S05A motor - not compatible with CL008-S05B For all CL009	Y		
CL008/15-S04B	CL008, CL008N & CL015	Main Assembly	Silica gasket (new)		Replace when usage is visible Includes height adjustment piece For CL008 from SN 13090130 on Not compatible with older models For all CL008N and for all CL015	Y		
CL008/9N-C-B	CL008, CL008N, CL009 & CL009N	Accessories & Consumables	Carbon brush (pair) -new	<a href="#">Picture / Photo</a>	Replace every 500 working hours, Replace every 250h if use is >5-6h/days CL008 from SN 1409.... on For CL008-S05B motor - not compatible with CL008-S05A For all CL008N For all CL009N	Y		
CL008/9-S06-110	CL008, CL008N, CL009 & CL009N	Main Assembly	Lid lock - 110V		From SN 1307... a new lid lock is used. The old and new model are interchangeable. Includes height adjustment piece	N		
CL008/9-S06-220	CL008, CL008N, CL009 & CL009N	Main Assembly	Lid lock - 220V		From SN 1307... a new lid lock is used. The old and new model are interchangeable. Includes height adjustment piece	N		
CL008/9-S07	CL008, CL008N, CL009 & CL009N	Main Assembly	Speed sensor			N		

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CL008/9-S19-110	CL008, CL008N, CL009 & CL009N	Main Assembly	Coil for AC filter 110V		Including capacitor	N		
CL008/9-S19-220	CL008, CL008N, CL009 & CL009N	Main Assembly	Coil for AC filter 220V		Including capacitor	N		
CF-T6.3AL250V	CL008, CL008N, CL009 & CL009N	Power	Fuse T6.3AL 250V 5x20mm (10 pcs)		At the back of CL008-110 and CL008N-110 At the back of CL009-220 and CL009N-220 On motor driver board of CL008-220	N		1:10
CF-T8.0AL250V	CL008, CL008N, CL009 & CL009N	Power	Fuse T8.0AL 250V 5x20mm (10 pcs)		On motor driver board of CL008-110 At the back of CL009-110 and CL009N-110	N		1:10
CL008/15-S10	CL008, CL008N, CL009, CL009N & CL015	Accessories & Consumables	Tube adaptors 5-7 ml (8 pcs)			N		
CL008/9/15-S09	CL008, CL008N, CL009, CL009N & CL015	Accessories & Consumables	Emergency opening bar			N		
CL008/9/15-S03	CL008, CL008N, CL009, CL009N & CL015	Main Assembly	Hinge		2 per instrument Includes height adjustment piece	N		
CL008/9/15-S12	CL008, CL008N, CL009, CL009N & CL015	Main Assembly	Motor mount (1 pcs)		6 per instrument Replace when motor has become unbalanced Recommended to change all 6 together	(Y)		
CL008/9/15-S13	CL008, CL008N, CL009, CL009N & CL015	Main Assembly	Fan			N		
CL008/9/15-S15	CL008, CL008N, CL009, CL009N & CL015	Main Assembly	Rubber foot/support (1 pcs)		4 per instrument	N		
CL008/9/15-S17	CL008, CL008N, CL009, CL009N & CL015	Main Assembly	Height adjustment for lid		Upon replacing the lidlock, the gasket or the hinges, it might be necessary to (slightly) adjust the height of the lid.	N		
CL008/9-S16-110	CL008, CL008N, CL009, CL009N & CL015	Main Assembly	Break resistor 110V			N		
CL008/9-S16-220	CL008, CL008N, CL009, CL009N & CL015	Main Assembly	Break resistor 220V			N		
CL008N-S01-110	CL008N	Main Assembly	Motor driver PC board - 110V		Only for CL008N-110 Not compatible with CL008	N		
CL008N-S01A-220	CL008N	Main Assembly	Motor driver PC board type A - 220V	<a href="#">Picture / Photo</a>	Only for CL008N, for SN 170200xx Not compatible with CL008	N		

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CL008N-S01B-220	CL008N	Main Assembly	Motor driver PC board type B - 220V	<a href="#">Picture / Photo</a>	Only for CL008N, from SN 1703xxxx on Not compatible with CL008	N		
CL008N-S02-110	CL008N	Main Assembly	DISPLAY - RPM/TIME PC board - 110V		Only for CL008N-110 Not compatible with CL008-110	N		
CL008N-S02-220	CL008N	Main Assembly	DISPLAY - RPM/TIME PC board -220V	<a href="#">Picture / Photo</a>	Only for CL008N-220 Not compatible with CL008-220	N		
CL008N-S18	CL008N	Main Assembly	Front label		Do NOT reattach the old label with (super)glue. Reattach with doublesided tape or use a new front label.	N		
CL008/15-S21	CL008N, CL009N, CL015	Accessories & Consumables	Tube adaptors 10-15mL (8pcs)			N		
CL008N/9N/15-S08	CL008N, CL009N, CL015	Main Assembly	Power switch		Not compatible with CL008 or CL009	N		
CL009-S10	CL009	Accessories & Consumables	Rubber for rotor (black)		Replace when usage is visible	Y		
CL009-S01-220	CL009	Main Assembly	Motor driver PC board (old) 220V		For SN untill 1204120 Not compatible with newer models	N		
CL009-S02A-220	CL009	Main Assembly	PC board for display (old, 2 pcs) 220V		For SN starting from 1012448 Not compatible with newer models	N		
CL009-S02B-110	CL009	Main Assembly	Display & Motor driver PC board (1 pc)			N		
CL009-S02B-220	CL009	Main Assembly	Display & Motor driver PC board (new,		For SN from 1102000 on Not compatible with older models	N		
CL009-S04A	CL009	Main Assembly	Rubber gasket for lid (old, red)		Replace when usage is visible Until SN 1308..... Not compatible with newer models Includes height adjustment piece	Y		
CL009-S04B	CL009	Main Assembly	Silica gasket for lid (new, red)		Replace when usage is visible From SN 13090130 on Not compatible with older models Includes height adjustment piece	Y		
CL009-S05-110	CL009	Main Assembly	Motor - 110V		12.000 rpm Only for CL009-110 Not compatible with CL009N-110	N		
CL009-S05-220	CL009	Main Assembly	Motor - 220V		12.000 rpm Only for CL009-220 Not compatible with CL009N-220	N		
CL009-S11	CL009	Main Assembly	Rotor			N		
CL009-S14	CL009	Main Assembly	Transformer for fan			N		
CL009-S18	CL009	Main Assembly	Front label		Do NOT reattach the old label with (super)glue. Reattach with doublesided tape or use a new front label.	N		

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CL009N-S02-110	CL009N	Main Assembly	Digital display & Motor driver PC board		Only for CL009N-110 Not compatible with CL009-110	N		
CL009N-S02-220	CL009N	Main Assembly	Digital display & Motor driver PC board	<a href="#">Picture / Photo</a>	Only for CL009N-220 Not compatible with CL009-220	N		
CL009N-S05-110	CL009N	Main Assembly	Motor - 110V		12.000 rpm Only for CL009N-110 Not compatible with CL009-110	N		
CL009N-S05-220	CL009N	Main Assembly	Motor - 220V	<a href="#">Picture / Photo</a>	12.000 rpm Only for CL009N-220 Not compatible with CL009-220	N		
CL009N-S18	CL009N	Main Assembly	Front label		Do NOT reattach the old label with (super)glue. Reattach with doublesided tape or use a new front label.	N		
CL010-C	CL010	Accessories & Consumables	Flat cover			N		
CL010-R	CL010	Accessories & Consumables	Test tube rack 50 places			(Y)		
CL010-S02	CL010	Main Assembly	Sensor			N		
CL010-S03-110	CL010	Main Assembly	Heater- 110V			N		
CL010-S03-220	CL010	Main Assembly	Heater- 220V	<a href="#">Picture / Photo</a>		N		
CL010-S04	CL010	Main Assembly	PID controller -	<a href="#">Picture / Photo</a>	Includes new front label	N		
CL010-S09	CL010	Main Assembly	Front label		After PCB or PID replacement, reattach with doublesided tape (NOT super glue) or use a new front label. New front label available on request with order of PCB or PID (mention SN)	N		
CF-CB-T10.0AL250V	CL010	Power	Circuit breaker T10.0AL250V		Used for CL010-110 from S/N 1904.... on Compatible with older instruments	N		1:10
CF-CB-T5.0AL250V	CL010	Power	Circuit breaker T5.0AL250V		For CL010-220	N		1:10
CF-CB-T4.0AL250V	CL010, CL011	Power	Circuit breaker T4.0AL250V		Used for CL010-220 until S/N 1902.... Use CF-CB-T5.0AL250V in case of replacement For CL011-110 and CL011-220	N		1:10
CL010/11/12-S01CABLE	CL010, CL011 & CL012	Main Assembly	Over-temperature connection cable for		When replacing PCB type A with PCB type B, C or D, and you want to add the over-temperature safety, it is necessary to order and install this cable as well	N		
CL010/11/12-S01D-110	CL010, CL011 & CL012	Main Assembly	PC Board - 110V (type D)		Type D with audible alarm for over-temperature. For CL011 Compatible with CL010 and CL012(N)	N		
CL010/11/12-S01D-220	CL010, CL011 & CL012	Main Assembly	PC Board - 220V (type D)	<a href="#">Picture / Photo</a>	Type D with audible alarm for over-temperature. For CL011 Compatible with CL010 and CL012(N)	N		

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CL010/11/12-S07	CL010, CL011 & CL012	Main Assembly	Solid state relay			N		
CL010/11/12-S10	CL010, CL011 & CL012	Main Assembly	EMI filter			N		
CL010/11/12-S01C-110	CL010, CL011 & CL014	Main Assembly	PC Board - 110V (type C)		Type C with over-temperature connection. For CL010 and CL012(N) Compatible with CL011	N		
CL010/11/12-S01C-220	CL010, CL011 & CL015	Main Assembly	PC Board - 220V (type C)	<a href="#">Picture / Photo</a>	Type C with over-temperature connection. For CL010 and CL012(N) Compatible with CL011	N		
CF-CB-T8.0AL250V	CL010, CL012 & CL012N	Power	Circuit breaker T8.0AL250V		Used for CL010-110 until S/N 1902.... Use CF-CB-T10.0AL250V in case of replacement For CL012(N)-220	N		1:10
CL011-S02	CL011	Main Assembly	Thermostat	<a href="#">Picture / Photo</a>		N		
CL011-S04-110	CL011	Main Assembly	Heater - 110V			N		
CL011-S04-220	CL011	Main Assembly	Heater - 220V	<a href="#">Picture / Photo</a>		N		
CL011-S05	CL011	Main Assembly	PID controller		Includes new front label	N		
CL011-S06	CL011	Main Assembly	Temperature sensor			N		
CL011-S09	CL011	Main Assembly	Front label		After PCB or PID replacement, reattach with doublesided tape (NOT super glue) or use a new front label. New front label available on request with order of PCB or PID (mention SN)	N		
CL011/12-S08	CL011 & CL012	Main Assembly	Rubber gasket (red)		Replace when usage is visible	Y		
CL012-S02	CL012 & CL012N	Main Assembly	Thermostat	<a href="#">Picture / Photo</a>		N		
CL012-S04-110	CL012 & CL012N	Main Assembly	Heater - 110V			N		
CL012-S04-220	CL012 & CL012N	Main Assembly	Heater - 220V	<a href="#">Picture / Photo</a>		N		
CL012-S05	CL012 & CL012N	Main Assembly	PID controller		Only for CL012 Not compatible with CL012N Includes new front label	N		
CL012-S06	CL012 & CL012N	Main Assembly	Temperature sensor			N		
CL012-S09	CL012	Main Assembly	Front label of CL012		After PID or PC board replacement Do NOT reattach the old label with (super)glue. Reattach with doublesided tape or use a new front label.	N		
CF-CB-T15.0AL250V	CL012	Power	Circuit breaker T15.0AL 250V		For CL012(N)-110	N		1:10
CL012N-S05	CL012N	Main Assembly	PID controller with timer function		Only for CL012N Not compatible with CL012	N		
CL012N-S09	CL012N	Main Assembly	Front label of CL012N		After PID or PC board replacement Do NOT reattach the old label with (super)glue. Reattach with doublesided tape or use a new front label.	N		
CL014-C1	CL014	Accessories & Consumables	Magnetic stirrer bar 25 x 8 mm			(Y)		

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CL014-C2	CL014	Accessories & Consumables	Magnetic stirrer bar 40 x 8 mm			(Y)		
CL015-S01	CL015	Main Assembly	PC Board	<a href="#">Picture / Photo</a>		N		
CL015-S05	CL015	Main Assembly	Motor	<a href="#">Picture / Photo</a>		N		
CL015-S06	CL015	Main Assembly	Lid lock		Includes height adjustment piece	N		
CL015-S18	CL015	Main Assembly	Front label		Do NOT reattach the old label with (super)glue. Reattach with doublesided tape or use a new front label.	N		
CL015-S20	CL015	Power	Power supply	<a href="#">Picture / Photo</a>		N		
CL016-S01	CL016	Power	Power adaptor		110-240VAC, 50/60 Hz	N		
CM001-S01	CM001	Accessories & Consumables	Halogen Lamp 6V-20W		Only for CM001 - Halogen microscope	Y		
CM001-S09	CM001	Accessories & Consumables	Bright field condenser		Standard part of CM001 - Halogen microscope	N		
CM001-S15	CM001	Main Assembly	Electronic boards		Only for CM001 - Halogen microscope	N		
CM001-S16	CM001	Main Assembly	Ceramic socket for lamp		Only for CM001 - Halogen microscope	N		
CM001-S17B	CM001	Main Assembly	Potentiometer (1 MΩ) - Type B		Only for CM001 - Halogen microscope	N		
CM001-S17C	CM001	Main Assembly	Potentiometer assembly (1 MΩ) - Type		Only for CM001 - Halogen microscope	N		
CM001-S18	CM001	Main Assembly	Transformer		Only for CM001 - Halogen microscope	N		
CM001-S02	CM001 & CM001N	Accessories & Consumables	EW 10x - 20mm			N		
CM001-S03	CM001 & CM001N	Accessories & Consumables	Infinitive semi plan objective 4x	<a href="#">Picture / Photo</a>		(Y)		
CM001-S04	CM001 & CM001N	Accessories & Consumables	Infinitive semi plan objective 10x	<a href="#">Picture / Photo</a>		(Y)		
CM001-S05	CM001 & CM001N	Accessories & Consumables	Infinitive semi plan objective 40x	<a href="#">Picture / Photo</a>		(Y)		
CM001-S06	CM001 & CM001N	Accessories & Consumables	Infinitive semi plan objective 100x	<a href="#">Picture / Photo</a>		(Y)		
CM001-S07	CM001 & CM001N	Accessories & Consumables	Nosepiece			N		
CM001-S08	CM001 & CM001N	Accessories & Consumables	Stage: 140x140mm			N		
CM001-S11A	CM001 & CM001N	Accessories & Consumables	Oil dark field condenser	<a href="#">Picture / Photo</a>	Optional accessory	N		
CM001-S12	CM001 & CM001N	Accessories & Consumables	Blue filter			N		
CM001-S13	CM001 & CM001N	Accessories & Consumables	Green filter			N		
CM001-S19	CM001 & CM001N	Accessories & Consumables	Handle for stage			N		
CM001N-S01A	CM001N	Accessories & Consumables	LED lamp - Type A	<a href="#">Picture / Photo</a>	Short line cable only for CM001N - LED microscope Until SN 103592	(Y)		
CM001N-S01B	CM001N	Accessories & Consumables	LED lamp - Type B	<a href="#">Picture / Photo</a>	Long line cable for CM001N - LED microscope Starting from SN 103593	(Y)		
CM001N-S09	CM001N	Accessories & Consumables	Bright field condenser		Standard part of CM001N - LED microscope	N		
CM001N-S15A	CM001N	Main Assembly	Electronic board - Type A	<a href="#">Picture / Photo</a>	Only for CM001N - LED microscope (from SN 103216 until SN 103592)	N		
CM001N-S15B	CM001N	Main Assembly	Electronic board - Type B	<a href="#">Picture / Photo</a>	Only for CM001N - LED microscope (from SN 103593 on)	N		
CM001N-S17A	CM001N	Main Assembly	Potentiometer (5kΩ) - Type A		Only for CM001N - LED microscope (from SN 103216 until SN 103592)	N		

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CM001N-S17B	CM001N	Main Assembly	Potentiometer assembly (2kΩ) - Type I	<a href="#">Picture / Photo</a>	Only for CM001N - LED microscope (from SN 103593 on)	N		
CM001N-S19	CM001N	Main Assembly	Dimming Circuit board	<a href="#">Picture / Photo</a>	Only for CM001N - LED microscope (starting from SN 103593)	N		
CM001N-S20	CM001N	Main Assembly	LED power module	<a href="#">Picture / Photo</a>	Only for CM001N - LED microscope (from SN 103216 until SN 103592)	N		
CM001N-S21	CM001N	Power	Power adapter		Only for CM001N - LED microscope (from SN 103593 on)	N		
CY003-S02A	CY003	Power	Power adaptor 18V (110-220V)			N		
CY003-S02B	CY003	Power	Power adaptor 18V 110V			N		
CY004-S10	CY004	Hydraulic system	Peristaltic pump with motor		Discontinued	N		1:10
CY004-S25	CY004	Hydraulic system	Pump tube		Replace every 6 months	Y	6M	1:5
CY004-S29A	CY004	Hydraulic system	Connector of waste tube (old)		For CYANStart until P/D 05/2014 Replace every 2-3 years	(Y)	2Y	1:5
CY004-S29B	CY004	Hydraulic system	Connector of waste tube (new)		From CYANStart P/D 05/2014 on Replace every 2-3 years	(Y)	2Y	1:5
CY004-S30	CY004	Hydraulic system	Connector of pump tube (1 unit)		2 per instrument Replace every 2-3 years	(Y)	2Y	1:3
CY004-S16	CY004	Main Assembly	Fan for power assembly (80x80)			N		1:10
CY004-S17	CY004	Main Assembly	Processor for main board			N		
CY004-S20	CY004	Main Assembly	Display (240*128)			N		1:10
CY004-S38	CY004	Main Assembly	Keyboard			N		1:10
CY004-S39	CY004	Main Assembly	Printerboard (w/o printer)		Only for old printers (CYANStart produced before 09/2014)	N		0
CY004-S01	CY004	Optic Assembly	Optical Assembly		correct alignment with motor, without filters, lamp, flowcell, heater, sensors, fan or board	N		
CY004-S11	CY004	Optic Assembly	Lamp (with frame, cable and connecto		Connection on power supply Replace every 2000 hours/ year	Y	1Y	1:2
CY004-S12	CY004	Optic Assembly	Motor of filter wheel			N		1:10
CY004-S15	CY004	Optic Assembly	Fan of optical assembly (50x50)			N		1:10
CY004-S33	CY004	Optic Assembly	Pre-amplifier board with cable		Replace every 2-3 years	(Y)	2Y	1:5
CY004-S34	CY004	Optic Assembly	Temperature sensor (optical)			N		1:10
CF-F2.0AL250V	CY004	Power	Fuse F2.0 AL250V 5x20mm (10 pcs)		2 per instrument	(Y)		1:10
CY004-S46	CY004	Power	Power plug connector (seat)			N		1:10
CY004/9-S43	CY004 & CY009	Accessories & Consumables	Dust cover			N		1:10
CY004/9-S27	CY004 & CY009	Hydraulic system	Metal aspiration tube		See hydraulic diagram	N		1:2
CY004/9-S44	CY004 & CY009	Hydraulic system	Waste bottle			(Y)		1:10
CY004/9-S21	CY004 & CY009	Main Assembly	Aspiration switch with cable	<a href="#">Picture / Photo</a>		N		1:10
CY004/9-S22	CY004 & CY009	Main Assembly	Aspiration button (black/blue)			N		1:10
CY004/9-S23	CY004 & CY009	Main Assembly	Printer with board			N		1:20



Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CY004/9-S50	CY004 & CY009	Main Assembly	Rubber foot/support		4 per instrument	N		1:10
CY009-S55	CY004 & CY009	Main Assembly	USB Socket with cable			N		1:10
CY004/9-S02	CY004 & CY009	Optic Assembly	Filter 340 nm			N		1:5
CY004/9-S03	CY004 & CY009	Optic Assembly	Filter 405 nm			N		1:10
CY004/9-S04	CY004 & CY009	Optic Assembly	Filter 492 nm			N		1:10
CY004/9-S05	CY004 & CY009	Optic Assembly	Filter 510 nm			N		1:10
CY004/9-S06	CY004 & CY009	Optic Assembly	Filter 546 nm			N		1:10
CY004/9-S07	CY004 & CY009	Optic Assembly	Filter 578 nm			N		1:10
CY004/9-S08	CY004 & CY009	Optic Assembly	Filter 620 nm			N		1:10
CY004/9-S09	CY004 & CY009	Optic Assembly	Flowcell			(Y)		1:5
CY004/9-S14	CY004 & CY009	Optic Assembly	Heater (Peltier) of flow cell			N		1:5
CY004/9-S35	CY004 & CY009	Optic Assembly	Home sensor of filterwheel			N		1:5
CY004/9-S45	CY004 & CY009	Optic Assembly	Fixing screw for lamp			N		1:10
CY004/9-S18	CY004 & CY009	Power	Power switch			N		1:10
CY004/9-S49	CY004 & CY009	Power	Power cord			N		1:10
CY006-MIBTR	CY006	Accessories & Consumables	CYANHemato Maintenance kit	<a href="#">Picture / Photo</a>		Y	2Y	1:2
CY006-S45	CY006	Accessories & Consumables	Tube adaptor 1, primary tubes			N		
CY006-S46	CY006	Accessories & Consumables	Tube adaptor 2, microtube			N		
CY006-S47	CY006	Accessories & Consumables	Tube adaptor 3, control vial			N		
CY006-S12	CY006	Cables	External Power Supply			N		
CY006-S49	CY006	Cables	Printer cable			N		
CY006-S50	CY006	Cables	Valves cable			N		
CY006-S51	CY006	Cables	LCD flat cable (white)			N		
CY006-S52	CY006	Cables	Backlight Touchscreen Start cable (gre			N		
CY006-S53	CY006	Cables	Preamplifier DIGIO cable (grey)			N		
CY006-S54	CY006	Cables	Preamplifier measure cable			N		
CY006-S55	CY006	Cables	Power cable (internal)			N		
CY006-S03	CY006	Dilutor Assembly	Micro piston			N		
CY006-S04	CY006	Dilutor Assembly	Macro piston			N		
CY006-S05	CY006	Dilutor Assembly	Lyse piston			N		
CY006-S37	CY006	Dilutor Assembly	Micro dilutor syringe	<a href="#">Picture / Photo</a>	Includes 3 x CY006-S07 O-rings	(Y)	3Y	1:5
CY006-S38	CY006	Dilutor Assembly	Macro dilutor syringe			N		1:5
CY006-S39	CY006	Dilutor Assembly	Lyse syringe			N		1:5
CY006-S13	CY006	Electronic Assembly	Amplifier board	<a href="#">Picture / Photo</a>		N		
CY006-S15	CY006	Electronic Assembly	Printer module		Without cable	N		
CY006-S31	CY006	Electronic Assembly	Touchscreen		A new touchscreen is used from SN 670538 on. When replacing an old touchscreen with a new it is necessary to replace the LCD support plate (CY006-S31S) as well.	N		1:5

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CY006-S31S	CY006	Electronic Assembly	Touchscreen support plate		When replacing an old touchscreen (SN until 670537), it is necessary to replace the LCD support plate (CY006-S31S) as well.	N		
CY006-S43	CY006	Electronic Assembly	DimmBoard	<a href="#">Picture / Photo</a>	When replacing a DimmPC (discontinued) with a DimmBoard, it is necessary to replace the mainboard as well	N		1:5
CY006-S44	CY006	Electronic Assembly	Main board (without) DimmBoard	<a href="#">Picture / Photo</a>	When replacing an old mainboard with DimmPC (discontinued) with a new mainboard, it is necessary to install a DimmBoard as well	N		
CY006-S60	CY006	Electronic Assembly	Micro SD card			N		
CY006-S11	CY006	Hydraulic system	Reagent tubing kit (tubes and caps)	<a href="#">Picture / Photo</a>	(includes 3xCY006-S16, 1xCY006-S23, 4xCY006-S21, 4xCY006-S26)	Y	2Y	1:5
CY006-S19	CY006	Hydraulic system	Tube to puffer reservoir and to pressu		2 tubes per instrument	(Y)	2Y	1:5
CY006-S20	CY006	Hydraulic system	Pump tube kit	<a href="#">Picture / Photo</a>		(Y)	2Y	1:5
CY006-S22	CY006	Hydraulic system	Internal hydraulic tubing (0,5 m)		You need 4 pcs on top of CY006-MIBTR to replace all tubes in the instrument To be cut to size	Y	2Y	1:5
CY006-S24	CY006	Hydraulic system	Big T connector (3,2 mm)		4 pcs per instrument	(Y)	3Y	1:5
CY006-S25	CY006	Hydraulic system	Small T connector (2,3 mm)		1 pcs per instrument	(Y)	3Y	1:5
CY006-S56	CY006	Hydraulic system	Steel reagent inlet connector (1 pcs)		4 pcs per instrument	(Y)	5Y	1:10
CY006-S59	CY006	Hydraulic system	Tube on needle head (1 cm)		Transfer or replace when replacing the needle or internal hydraulic sampling tube	(Y)		1:5
CY006-S17	CY006	Motor Assembly	Belt needle			(Y)	2Y	1:5
CY006-S18	CY006	Motor Assembly	Belt sample rotor			(Y)	2Y	1:5
CY006-S30	CY006	Motor Assembly	Opto cable (grey)			N		1:5
CY006-S32	CY006	Motor Assembly	Optosensor (black)			N		1:5
CY006-S40	CY006	Motor Assembly	Stepper motor			N		
CY006-S09	CY006	Pressure and Vacuum	Pump		Including electronic board and cable	N		
CY006-S28	CY006	Pressure and Vacuum	Electrovalve 2/2	<a href="#">Picture / Photo</a>	For V1, V6, V7	N		1:5
CY006-S29	CY006	Pressure and Vacuum	Electrovalve 3/2	<a href="#">Picture / Photo</a>	For V2, V3, V4, V5, V8, V9, V10	N		1:5
CY006-S41	CY006	Pressure and Vacuum	Puffer reservoir			N		
CY006-S42	CY006	Pressure and Vacuum	Puffer reservoir holder (1 pcs)		2 per instrument	N		
CY006-S57	CY006	Pressure and Vacuum	Valve block 1-5		Including electrovalve board and valves V1-V5	N		
CY006-S58	CY006	Pressure and Vacuum	Valve block 6-10		Including electrovalve board and valves V6-V10	N		
CY006-S02	CY006	Sampling & Measuring Assem	Sealing ring for microdilutor (1 pcs)	<a href="#">Picture / Photo</a>		(Y)	2Y	1:3
CY006-S10	CY006	Sampling & Measuring Assem	HGB measuring head	<a href="#">Picture / Photo</a>		N		1:5
CY006-S14	CY006	Sampling & Measuring Assem	Aperture 70 µm	<a href="#">Picture / Photo</a>	Includes CY006-S06 O-ring Replace every 3-5 years	(Y)	3Y	1:5
CY006-S27	CY006	Sampling & Measuring Assem	Needle		Replace every 3-5 years	(Y)	3Y	1:5

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CY006-S35	CY006	Sampling & Measuring Assembl	Reference electrode	<a href="#">Picture / Photo</a>	Includes CY006-S06 O-ring Replace every 3-5 years	(Y)	3Y	1:5
CY006-S36	CY006	Sampling & Measuring Assembl	Counting Chamber		Replace every 3-5 years	(Y)	3Y	1:5
CY006-S48	CY006	Sampling & Measuring Assembl	U-shaped metal fixing for aperture			N		
CY005/8-C21	CY008	Accessories & Consumables	0,8 ml cups - Sample /Pre-dilution (20l		For predilution	Y	User dependent	User dependent
CY008-C01	CY008	Accessories & Consumables	Drying pad		Replace every 6 months	Y	2Y	1:4
CY008-C02	CY008	Accessories & Consumables	Rack for primary tubes (universal) + 3r			N		
CY008-C03	CY008	Accessories & Consumables	Reaction racks cover			N		
CY008-C04	CY008	Accessories & Consumables	Rack for predilution cups (0,8 ml)		Optional feature	N		2:5
CY008-C06	CY008	Accessories & Consumables	Bottle assembly 2L with drilled cap			(Y)		1:5
CY008-C08	CY008	Accessories & Consumables	Reagent tray cover			N		
CY008-C09	CY008	Accessories & Consumables	Reagent tray (without cover)		Sold as compatible pair Follow the procedure to check the alignment when using different reagent trays	N		
CY008-C10	CY008	Accessories & Consumables	Tubing kit peristaltic pump (2 pcs)		Replace every 6 months	Y	6M	1:1
CY008-C11	CY008	Accessories & Consumables	Tubing kit - complete		Replace every year	Y	1Y	1:2
CY008-C12	CY008	Accessories & Consumables	Fitting kit red (1 pcs)		Replace every 2-3 years	(Y)	2Y	1:5
CY008-C13	CY008	Accessories & Consumables	Fitting kit blue (1 pcs)		Replace every 2-3 years	(Y)	2Y	1:5
CY008-C14	CY008	Accessories & Consumables	Fitting kit yellow (1 pcs)		Replace every 2-3 years	(Y)	2Y	1:5
CY008-C15	CY008	Accessories & Consumables	Fitting kit green (1 pcs)		Replace every 2-3 years	(Y)	2Y	1:5
CY008-C17	CY008	Accessories & Consumables	Inlet tubes kit (4 pcs)		Replace every year	(Y)	1Y	2:5
CY008-C20	CY008	Accessories & Consumables	Halogen lamp (6 V - 10 W)		Replace every 2000 hours/ year	Y	1Y	2:5
CY008-C22	CY008	Accessories & Consumables	1 ml cups - Calibrator & Control (100C		For calibrator and controls	Y		2:5
CY008-C25	CY008	Accessories & Consumables	Reagent containers 35 ml + caps (12 p			Y		2:5
CY008-C26	CY008	Accessories & Consumables	Reagent container 6 mL + caps (12 pcs			Y		2:5
CY008-C27	CY008	Accessories & Consumables	Reaction cuvettes (6 sectors - 120 pcs)		Replace every 3-6 months	Y		2:5
CY008-C28	CY008	Accessories & Consumables	3 ml cups - Samples in primary tube ra		For samples in primary tube rack	Y		2:5
CY008-C29	CY008	Accessories & Consumables	Kit for T° adjustment & calibration			N		
CY008-C31	CY008	Accessories & Consumables	Probe cleaning wire		Length 170 mm	N		1:5
CY008-C32	CY008	Accessories & Consumables	Tool for peristaltic pump adjustment			N		
CY008-S30	CY008	Arm Assembly	Sampling probe (internal cannula)			N		
CY008-S32	CY008	Arm Assembly	Level sensor and pre-heater assembly			N		
CY008-S33	CY008	Arm Assembly	Sampling probe assembly (internal anc			N		
CY008-S34	CY008	Arm Assembly	Arm assembly (whole)			N		
CY008-S35	CY008	Arm Assembly	Belt - arm vertical movement			N		
CY008-S36	CY008	Arm Assembly	Belt - arm horizontal movement			N		
CY008-S37	CY008	Arm Assembly	Home sensor assembly (Arm H/V, Wa			N		
CY008-S38	CY008	Arm Assembly	Pre- heater fitting kit (1 pcs)			N		
CY008-S80	CY008	Arm Assembly	Special key for probe tightening: dyna			N		
CY008-S31	CY008	Arm Assembly, Hydraulic Asse	Stepper motor assembly (Arm H/V, O			N		

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CY005/8-S018	CY008	Cables	External Serial Cable (m/f 9 pin, 5m)		For the most stable connexion	N		
CY005/8-S019	CY008	Cables	Power cord (external)			N		
CY005/8-S93	CY008	Cables	RS232 to USB converter cable with sof		Preferally use CY005/8-S018 to make the connection with the PC	N		
CY008-S81	CY008	Cables	Flat cable 10 PINCY008-S14 to CY008-S		See service manual chapter 3 for schemes of the cables	N		
CY008-S82	CY008	Cables	Flat cable 16 PINCY008-S15 to CY008-S		See service manual chapter 3 for schemes of the cables	N		
CY008-S83	CY008	Cables	Flat cable 20 PINCY008-S15 to CY008-S		See service manual chapter 3 for schemes of the cables	N		
CY008-S84	CY008	Cables	Flat cable 10 PINCY008-S14 to CY008-S		See service manual chapter 3 for schemes of the cables	N		
CY008-S85	CY008	Cables	Flat cable 24 PINCY008-S14 to CY008-S		See service manual chapter 3 for schemes of the cables	N		
CY008-S86	CY008	Cables	Flat cable 26 PINCY008-S14 to CY008-S		See service manual chapter 3 for schemes of the cables	N		
CY008-S87	CY008	Cables	Flat cable 34 PINCY008-S14 to CY008-S		See service manual chapter 3 for schemes of the cables	N		
CY008-S88	CY008	Cables	Internal Serial cable		See service manual chapter 3 for schemes of the cables	N		
CY008-S89	CY008	Cables	Flat cable 10 PINCY008-S17 to CY008-S		See service manual chapter 3 for schemes of the cables	N		
CY008-S90	CY008	Cables	Bottles connecting cable (black)			N		
CY008-S91	CY008	Cables	Flat cable 10 PINCY008-S13 to CY008-S		See service manual chapter 3 for schemes of the cables	N		
CY005/8-S001	CY008	Dilutor Assembly	Diluter assembly			N		
CY005/8-S002	CY008	Dilutor Assembly	Diluter fitting (1 pcs)		3 per instrument	(Y)	2Y	1:5
CY005/8-S044	CY008	Dilutor Assembly	Home sensor assembly Pro: A+D+P; Ex			N		
CY008-S39	CY008	Dilutor Assembly	Electrovalve (EV1 for diluter)			N		1:5
CY008-S40	CY008	Dilutor Assembly	Micro-pump assembly (µP1 for diluter			N		1:5
CY008-S43	CY008	Dilutor Assembly	Tubing kit for diluter (3 PCS)		Or part of CY008-C11 tubing kit	(Y)		1:5
CY008-S45	CY008	Dilutor Assembly	EV1 plastic fitting (1 pcs)			(Y)	2Y	1:5
CY008-S10	CY008	Electronic Assembly	Pre-amplifier board (sample channel)			N		
CY008-S11	CY008	Electronic Assembly	Pre-amplifier board (reference channe			N		
CY008-S12	CY008	Electronic Assembly	Plate interface board			N		1:5
CY008-S13	CY008	Electronic Assembly	Sampling arm interface board			N		1:5
CY008-S14	CY008	Electronic Assembly	Analytical control board			N		1:5
CY008-S15	CY008	Electronic Assembly	Stepper motors driver board			N		1:5
CY008-S16	CY008	Electronic Assembly	Photometer lamp board			N		
CY008-S17	CY008	Electronic Assembly	Electrovalve controller board			N		1:5

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CY008-S18	CY008	Electronic Assembly	Control panel board			N		
CY005/8-S005	CY008	Hydraulic system	Fitting nipples pump (1 pcs)		2 per instrument Replace every 2-3 years	(Y)		1:5
CY008-S29	CY008	Hydraulic system	Washing well assembly			N		
CY008-S46	CY008	Hydraulic system	Washing cannula assembly "1 & 2"			(Y)	2Y	1:5
CY008-S47	CY008	Hydraulic system	Washing cannula assembly "3"			(Y)	2Y	1:5
CY008-S48	CY008	Hydraulic system	Washing cannula assembly "4"			(Y)	2Y	1:5
CY008-S49	CY008	Hydraulic system	Washing can. assembly "5" w/o drying			(Y)	2Y	1:5
CY008-S51	CY008	Hydraulic system	Washing station probes assembly (whr			N		
CY008-S53	CY008	Hydraulic system	Spring for cannula (5 pcs)			N		
CY008-S54	CY008	Hydraulic system	Micro-pump assembly (µP 2 - µP 5)			N		1:5
CY008-S55	CY008	Hydraulic system	Micro-air pump assembly (µP 6)			N		1:5
CY008-S56	CY008	Hydraulic system	Peristaltic pump assembly			N		
CY008-S58	CY008	Hydraulic system	Solenoid valve - 2 way (EV2, EV3, EV4)			N		1:5
CY008-S59	CY008	Hydraulic system	Solenoid valve - 3 way (EV5, EV8)			N		1:5
CY008-S60	CY008	Hydraulic system	Peristaltic pump rotor			N		
CY008-S75	CY008	Hydraulic system	Inlet/Outlet fitting for WS-I, -II, -III and		4 per instrument or part of CY008-C17 inlet tubes kit Replace every 2-3 years	(Y)	2Y	1:5
CY008-S76	CY008	Hydraulic system	Bottle level sensor assembly			N		
CY008-S77	CY008	Hydraulic system	Bottle 2L			(Y)		1:4
CY008-S94	CY008	Hydraulic system	Filter after washing well			(Y)	2Y	1:5
CY008-S02	CY008	Main Assembly	Fan assembly (optic)			N		
CY008-S03	CY008	Main Assembly	Fan assembly (basement and rear side			N		
CY008-S05	CY008	Main Assembly	Cover lock switch			N		
CY008-S07	CY008	Main Assembly	Air spring for cover			N		
CY008-S63	CY008	Optic Assembly	Photometer assembly (2 optic channel			N		
CY008-S64	CY008	Optic Assembly	Interferential filters kit "Stablifite"			N		1:5
CY008-S65	CY008	Optic Assembly	Home sensor assembly (Optic + Plates			N		1:5
CY008-S66	CY008	Plates Assemblies	Upper pulley reaction plate			N		
CY008-S67	CY008	Plates Assemblies	Middle pulley sample plate			N		
CY008-S68	CY008	Plates Assemblies	Lower pulley reagent plate			N		
CY008-S69	CY008	Plates Assemblies	Ambient temperature sensor assembly			N		
CY008-S70	CY008	Plates Assemblies	Reaction plate heater belt			N		
CY008-S71	CY008	Plates Assemblies	Motor assy with pulley clamp (all plate			N		
CY008-S72	CY008	Plates Assemblies	Temperature sensor assembly			N		
CY008-S73	CY008	Plates Assemblies	Belt for plate (reaction, sample and re:			N		1:5
CY008-S04	CY008	Power	Main power socket			N		
CY008-S06A	CY008	Power	Power supply (old)		For SN 001 till SN 018. Not interchangeable with CY008-S06B	N		

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CY008-S06B	CY008	Power	Power supply (new)		From SN 019 on a new power supply is used. Not interchangeable with CY008-S06A	N		
CY008-S09	CY008	Power	Fuses T6.3AH250V 5x20mm (10 pcs)		2 per instrument	N		1:5
CY009-MIBTR	CY009	Accessories & Consumables	CYANSmart Maintenance kit	<a href="#">Picture / Photo</a>		Y	1Y	1:2
CY009-S10	CY009	Hydraulic system	Peristaltic pump rotor with motor			N		
CY009-S29	CY009	Hydraulic system	Connector of waste tube		See hydraulic diagram	(Y)	2Y	1:5
CY009-S30	CY009	Hydraulic system	Connector of pump tube (2 pcs)		See hydraulic diagram	(Y)	2Y	1:4
CY009-S54	CY009	Hydraulic system	Rubber of peristaltic pump (2 pcs)	<a href="#">Picture / Photo</a>		(Y)	2Y	1:5
CY009-S13	CY009	Main Assembly	Main Control Board			N		1:10
CY009-S16B	CY009	Main Assembly	Fans for optical and power assembly		From SN BS6KMxxx; P/D 2019/06 on Compatible with older instruments Connect on mainboard on connector B3	N		1:10
CY009-S19	CY009	Main Assembly	Cable of touch screen			N		1:10
CY009-S20A	CY009	Main Assembly	Touch screen (Type A)		For SN BS6GHxxx Software version 1.0.11B & 1.0.11 Includes CY009-S19 touchscreen cable	N		1:10
CY009-S20B	CY009	Main Assembly	Touch screen (Type B)		For all other SN Software version 1.1.1B or higher Includes CY009-S19 touchscreen cable	N		1:10
CY009-S34	CY009	Main Assembly	Temperature sensor (optical & incubat			N		1:10
CY009-S42	CY009	Main Assembly	Case			N		
CY009-S51	CY009	Main Assembly	Incubating block/ heating block			N		1:10
CY009-S52	CY009	Main Assembly	Heater element for incubator with ten			N		1:10
CY009-S56	CY009	Main Assembly	USB cable		For connection with PC software	N		1:10
CY009-S01	CY009	Optic Assembly	Optical Assembly		Correct alignment with filter wheel motor Does not include filters, lamp, flowcell, heater, sensors, fan or amplifier board	N		
CY009-S11	CY009	Optic Assembly	Lamp (with frame, cable and connecto	<a href="#">Picture / Photo</a>	Replace every 2000 hours/year Connection directly on main board (6V 10W)	Y	1Y	1:2
CY009-S12	CY009	Optic Assembly	Motor of filter wheel			N		1:10
CY009-S33	CY009	Optic Assembly	Pre-amplifier board with cable			Y	2Y	1:2
CY009-S24B	CY009	Power	Power supply (type B)	<a href="#">Picture / Photo</a>	From SN BS6Ilxxx; P/D 2017/03/01 on Compatible with older instruments Cables and support included	N		1:10
CY009-S46	CY009	Power	Power plug connector (seat)			N		1:10
CF-T2.0AL250V	CY009, CL001, CL015	Power	Fuse T2,0AL250V 5x20mm (10 pcs)		2 per instrument, at the backside (CY009) On mainboard (CL001) At the backside (CL015)	(Y)		1:10
CY010-S001	CY010	Accessories & Consumables	CYANCyte 560 Maintenance kit	<a href="#">Picture / Photo</a>		Y	1Y	1:5
CY010-S003	CY010	Accessories & Consumables	Needle cleaner wire		Or part of specific toolkit CY010-S005	N		1:5

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CY010-S005	CY010	Accessories & Consumables	CYANCyte 560 Specific toolkit	<a href="#">Picture / Photo</a>	Contains sample collar tool, UNF tools (2 types), UNF-UNF connector (2), tube connectors (2 types), needle cleaning wire (CY010-S003)	N		1:5
CY010-S007	CY010	Accessories & Consumables	Wifi Dongle		Wireless 11n 150 Mbps USB Adapter	N		
CY010-S009	CY010	Accessories & Consumables	Cleaning vials (3 pcs)		For needle cleaning process	Y		1:1
CY010-S078	CY010	Accessories & Consumables	Fuse F3.5A125V SMD (10 pcs)		On source board (S077)	(Y)		1:5
CY010-S089	CY010	Accessories & Consumables	External power supply		Power cord not included	N		
CY010-S046	CY010	Electronic Assembly	Pressure sensor board			N		
CY010-S070	CY010	Electronic Assembly	LILIOM board	<a href="#">Picture / Photo</a>	Battery is part of maintenance kit (S001) Or part of diplay panel module (S085)	N		
CY010-S071	CY010	Electronic Assembly	Speaker			N		
CY010-S072	CY010	Electronic Assembly	COLIBRI board	<a href="#">Picture / Photo</a>		N		
CY010-S073	CY010	Electronic Assembly	HVDCDC board	<a href="#">Picture / Photo</a>	High voltage generator board	N		
CY010-S074	CY010	Electronic Assembly	FCTRL board	<a href="#">Picture / Photo</a>	Fluidics control board	N		
CY010-S075	CY010	Electronic Assembly	Start button board	<a href="#">Picture / Photo</a>		N		
CY010-S076	CY010	Electronic Assembly	COMCON board			N		
CY010-S077	CY010	Electronic Assembly	Source board		Including cable to Start button board	N		
CY010-S079	CY010	Electronic Assembly	LVDS adapter board	<a href="#">Picture / Photo</a>	From display panel to LILIOM board Or part of diplay panel module (S085)	N		
CY010-S080	CY010	Electronic Assembly	IMP-AMP board	<a href="#">Picture / Photo</a>	Impedance amplifier control board	N		
CY010-S081	CY010	Electronic Assembly	D5ACQ board	<a href="#">Picture / Photo</a>	Data processing board (impedance & optical)	N		
CY010-S082	CY010	Electronic Assembly	Front panel camera			N		
CY010-S083	CY010	Electronic Assembly	Touch controller	<a href="#">Picture / Photo</a>		N		
CY010-S085	CY010	Electronics	Display panel module		Incl. LCD display, LVDS adapter (S079), flat display cable, LILIOM board (S070) and metal carrier plate	N		
CY010-S011	CY010	Hydraulic system	Sample vial holder		Including o-ring and screw Or part of sample rotor module (S012)	N		
CY010-S012	CY010	Hydraulic system	Sample rotor module		Including motor, rotor disk, sensor set (S014), vial holder (S011), micro-switch (S020), blood sample sensor (S023) Without FCTRL board	N		
CY010-S014	CY010	Hydraulic system	Sample rotor latch sensor set		Or part of sample rotor module (S012)	N		
CY010-S015	CY010	Hydraulic system	Ceramic valve module	<a href="#">Picture / Photo</a>	Including shear valve, needle (S016), washing head (S017), leak chamber (S018), blood sensor (S021), all connected tubes and SCTRL board	N		
CY010-S016	CY010	Hydraulic system	Piercing needle		With UNF and rubber o-ring Or part of yearly maintenance kit CY010-S001 or ceramic valve module (S015)	Y		1:5
CY010-S017	CY010	Hydraulic system	Needle washing head	<a href="#">Picture / Photo</a>	Or part of ceramic valve module (S015)	N		
CY010-S018	CY010	Hydraulic system	Needle leak collection tray		Or part of ceramic valve module (S015)	N		

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CY010-S019	CY010	Hydraulic system	Cable of needle washing head and lead			N		
CY010-S020	CY010	Hydraulic system	Piercing micro-switch	<a href="#">Picture / Photo</a>	Or part of sample rotor module (S012)	N		
CY010-S023	CY010	Hydraulic system	Blood sample sensor	<a href="#">Picture / Photo</a>	Or part of sample rotor module (S012)	N		
CY010-S024	CY010	Hydraulic system	Sampling system tube set	<a href="#">Picture / Photo</a>	Including all tubes related to the sampling system, with connectors	(Y)		1:5
CY010-S037	CY010	Hydraulic system	Micro injector		Including motor and syringe	N		
CY010-S050	CY010	Hydraulic system	External reagent tube set			Y		1:5
CY010-S051	CY010	Hydraulic system	Reagent connector module		Plate with 5 metal nozzles	N		
CY010-S055	CY010	Hydraulic system	Injector module	<a href="#">Picture / Photo</a>	Including motor, syringe and ICTRL board Without IFS or valves For FI, BI, WI, SI, LI, RI Neutral address	N		1:5
CY010-S056	CY010	Hydraulic system	Injector fluid sensor (IFS)		For FI, BI, WI, LI	N		
CY010-S057	CY010	Hydraulic system	RI injector fluid sensor (IFS)		Including cable and support	N		
CY010-S058	CY010	Hydraulic system	Valve (3/2)	<a href="#">Picture / Photo</a>	Including cable, connector and sealing element For valves E, T, Q, F, K, B, A, W, N, I, L, R, C	N		1:5
CY010-S059	CY010	Hydraulic system	Valve (2/2)	<a href="#">Picture / Photo</a>	Incl. cable, connector and sealing element For valve S	N		
CY010-S060	CY010	Hydraulic system	Valve manifold			(Y)		1:5
CY010-S062	CY010	Hydraulic system	Peristaltic pump tubes (2pcs)		Or part of yearly maintenance kit CY010-S001	(Y)	1Y	1:10
CY010-S063	CY010	Hydraulic system	Peristaltic pump			N		
CY010-S065	CY010	Hydraulic system	Internal reagent inlet tube set	<a href="#">Picture / Photo</a>	All tubes related to internal reagent inlet, with connectors	(Y)		1:5
CY010-S086	CY010	Power	Front panel cable set	<a href="#">Picture / Photo</a>		N		
CY010-S087	CY010	Power	LIN cable set	<a href="#">Picture / Photo</a>		N		
CY010-S088	CY010	Power	Internal power cable		From power socket to SOURCE board	N		
CY010-S030	CY010	RBC Channel	Impedance module		Including RBC aperture module (S031), micro injector (S037), RBC mixer (S034), mix chamber (S035), RCTRL board, valves, all connected tubes	N		
CY010-S031	CY010	RBC Channel	RBC aperture module (70 µm)	<a href="#">Picture / Photo</a>	Including aperture, duckbill valve, electrodes or part of impedance module (S030)	N		1:5
CY010-S033	CY010	RBC Channel	Electronic shield for RBC module			N		
CY010-S034	CY010	RBC Channel	RBC sample mixer with motor		Or part of Impedance module (S030)	N		
CY010-S035	CY010	RBC Channel	Mix and waste collection chamber		Or part of Impedance module (S030)	N		
CY010-S036	CY010	RBC Channel	RBC system tube set		All tubes related to the RBC channel, with connectors	(Y)		1:5
CY010-S040	CY010	WBC channel	Optical module	<a href="#">Picture / Photo</a>	Including DHGB and OPT AMP boards and all connected tubes	N		
CY010-S041	CY010	WBC channel	Meander temperature module	<a href="#">Picture / Photo</a>	Includng peltier, temperature sensor (S042) and radiator element	N		



Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CY010-S042	CY010	WBC channel	Meander temperature sensor		Or part of meander temperature module (S041)	N		
CY010-S043	CY010	WBC channel	WBC system tube set		All tubes related to WBC channel, with connectors	(Y)		1:5
CY010-S045	CY010	WBC channel	Pressure sensor			N		
CY011-S03	CY011	Accessories & Consumables	USB cable			N		
CY011-S05	CY011	Accessories & Consumables	Strip loader			(Y)		1:5
TR	CY011	Accessories & Consumables	CYANStrip Mini repair kit			N		
CY011-S09	CY011	Power	Power adaptor			N		1:5
CY012-S02	CY012	Power	Power adaptor 15V (110-220V)			N		
CY014-MIBTR	CY014	Accessories & Consumables	CYANVision Maintenance kit	<a href="#">Picture / Photo</a>		Y	1Y	1:2
CY014-S30	CY014	Electronics	BCAM (incl. Mainboard & ON/OFF Swi			N		
CY014-S31	CY014	Electronics	Optical system (excl. Mainboard, inclu			N		
CY014-S32	CY014	Electronics	Peltier element			N		1:10
CY014-S33	CY014	Electronics	Mainboard			N		
CY014-S34	CY014	Electronics	Main CPU battery			(Y)		1:5
CY014-S07	CY014	Hydraulic system	Peristaltic Pump (Incl. cassette & cable)			N		1:10
CY014-S11	CY014	Hydraulic system	Waste bottle			N		1:10
CY014-S20-01	CY014	Main Assembly	Printer			N		1:10
CY014-S20-02	CY014	Main Assembly	Printer data cable			N		
CY014-S20-03	CY014	Main Assembly	Printer power cable			N		
CY014-S21	CY014	Main Assembly	HDMI extension cable			N		
CY014-S22	CY014	Main Assembly	Push button (incl.cable)			N		
CY014-S36	CY014	Main Assembly	ON/OFF Button			N		
CY014-S01	CY014	Optic Assembly	Flowcell assembly			N		1:5
CY014-S03	CY014	Optic Assembly	Flowcell outlet connector			N		
CY014-S40	CY014	Power	Power supply (Incl. EU Power cord)			N		
CY015-C01	CY015	Accessories & Consumables	Fuse F2.5AL250V 5x20mm (5 pcs)			(Y)		1:5
CY015-C02	CY015	Accessories & Consumables	Halogen lamp (6V, 10W)	<a href="#">Picture / Photo</a>	Replace every 1000 working hours	Y	1000 hours	1:2
CY015-C04	CY015	Accessories & Consumables	2mL reagent cups 10,8x40mm (30 pcs)	<a href="#">Picture / Photo</a>		Y	User dependen	2:5
CY015-C05	CY015	Accessories & Consumables	1mL sample/reagent cups (1000 pcs)	<a href="#">Picture / Photo</a>		Y	User dependen	2:5
CY015-C06	CY015	Accessories & Consumables	15mL reagent bottle (30 pcs)	<a href="#">Picture / Photo</a>		Y	User dependen	2:5
CY015-C07	CY015	Accessories & Consumables	35mL reagent bottle (30 pcs)	<a href="#">Picture / Photo</a>		Y	User dependen	2:5
CY015-C08	CY015	Accessories & Consumables	Screw caps for reagents (30 pcs)			Y	User dependen	2:5
CY015-C09	CY015	Accessories & Consumables	Reaction segments of 24 cuvettes (10C	<a href="#">Picture / Photo</a>		Y	User dependen	1:2
CY015-S01	CY015	Accessories & Consumables	Adapter for 15mL reagent bottles (15	<a href="#">Picture / Photo</a>		N		1:10
CY015-S02	CY015	Accessories & Consumables	Adapter for 2mL or 1mL reagent cups	<a href="#">Picture / Photo</a>		N		1:10
CY015-S03	CY015	Accessories & Consumables	Adapter for 1mL sample cups (10 pcs)	<a href="#">Picture / Photo</a>		N		1:10
CY015-S04	CY015	Accessories & Consumables	2L container			N		1:10
CY015-S05	CY015	Accessories & Consumables	External tubing and level sensor assem			N		
CY015-S07	CY015	Accessories & Consumables	Sample tray (position 1-10)			N		

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CY015-S13	CY015	Accessories & Consumables	Peristaltic pump cassette	<a href="#">Picture / Photo</a>	2 per instrument Replace every 100 working hours	Y	100 hours	2:1
CY015-S33	CY015	Accessories & Consumables	Reagent tray (20 positions)			N		
CY015-S24	CY015	Arm assembly	Level sensor board			N		
CY015-S25	CY015	Arm assembly	Sampling needle holder with knob and			N		1:10
CY015-S26	CY015	Arm assembly	Level sensor contact			N		
CY015-S29	CY015	Arm assembly	Sampling arm belt			N		1:5
CY015-S71	CY015	Arm assembly	Motor for needle vertical			N		
CY015-S85	CY015	Arm assembly	Needle horizontal home sensor			N		
CY015-S86	CY015	Arm assembly	Needle vertical home sensor			N		
CY015-S87	CY015	Arm assembly	Level sensor cable			N		
CY015-S64	CY015	Electronic system	CPU board (on request, SN specific)			N		
CY015-S65	CY015	Electronic system	Power board			N		1:10
CY015-S66	CY015	Electronic system	Power supply			N		1:10
CY015-S88	CY015	Electronic system	Internal USB serial cable			N		
CY015-S89	CY015	Electronic system	Interboard cable			N		
CY015-S91	CY015	Electronic system	USB converter board			N		
CY015-S92	CY015	Electronic system	Bottom fan			N		
CY015-S11	CY015	Hydraulic system	Hydraulic panel assembly			N		
CY015-S12	CY015	Hydraulic system	Peristaltic pump motor and cassette			N		1:10
CY015-S15	CY015	Hydraulic system	Diluter assembly			N		
CY015-S16	CY015	Hydraulic system	Diluter gasket		Replace every year	Y	1Y	1:2
CY015-S17	CY015	Hydraulic system	Tubing diluter to needle + 1 fitting			(Y)		1:5
CY015-S18	CY015	Hydraulic system	3-way valve			N		
CY015-S19	CY015	Hydraulic system	Tubing valve to diluter + 2 fittings + sh			(Y)		1:5
CY015-S23	CY015	Hydraulic system	Sampling needle		Replace every year or in case of bending	Y	1Y	1:2
CY015-S73	CY015	Hydraulic system	Diluter motor			N		
CY015-S79	CY015	Hydraulic system	Tubing washing well to waste			(Y)		1:5
CY015-S41	CY015	Optic Assembly	Pre-amplifier board			N		
CY015-S42	CY015	Optic Assembly	Pre-amplifier cable			N		
CY015-S46	CY015	Optic Assembly	Filter wheel assembly (incl. filter set)			N		
CY015-S47	CY015	Optic Assembly	Filter wheel (without filters)			N		
CY015-S49	CY015	Optic Assembly	Complete optical assembly (incl. filter)			N		
CY015-S60	CY015	Optic Assembly	Lamp fastening knob			N		
CY015-S72	CY015	Optic Assembly	Filter wheel motor			N		
CY015-S74	CY015	Optic Assembly	Optical home sensor			N		
CY015-S84	CY015	Optic Assembly	Optical cable			N		
CY015-S36	CY015	Plates assembly	Complete cooling plate assembly (plas			N		
CY015-S37	CY015	Plates assembly	Peltier element for cooling plate			N		1:10
CY015-S38	CY015	Plates assembly	Fan for cooling plate 24V			N		

Price list data (dynamic)

Item Code	Analyzer	System	Description	See picture	Comments	Consumable	Replacement interval	Suggested stock (for analyzers only)
CY015-S43	CY015	Plates assembly	Cuvette heating resistance			N		1:10
CY015-S44	CY015	Plates assembly	Cuvette temperature sensor			N		1:10
CY015-S45	CY015	Plates assembly	Cuvette plate belt			N		
CY015-S70	CY015	Plates assembly	Motor for needle horizontal and optical			N		
CF-T4.0AL250V	CY015	Power	Fuse T4.0AL250V 5x20mm (10 pcs)		2 per instrument	(Y)		1:10
CY006-S08	CY006	Hydraulic system	Washing head		Also included in CY006-MIBTR	Y	1Y	1:5